

DISTRIBUTED BY



SOUTHEASTERN
BIOMEDICAL ASSOCIATES INC

3546 Hickory Blvd. Hudson, NC 28638


888.310.7322 www.sebiomedical.com

Email: info@sebiomedical.com

AARON® 950™

Bovie®

ELECTROSURGERY FOR THE OFFICE AND SURGERY CENTER



VERSATILITY

High Frequency Desiccator with 60 watts of Cutting capabilities

UNIQUE

Two machines in one convenient package

SIMPLICITY

Six (6) Doctor presets for fast, easy setup

ACCURACY

Adjustable in 1/10th watt increments in the first 10 watts, for the most accurate output available

COST EFFECTIVE

No grounding pad required when working in the desiccation mode.

RELIABILITY

The Aaron 950 is backed by Bovie's two (2) year manufacturer's limited warranty.

ELECTROSURGICAL GENERATOR





AARON® 950™ INCLUDES

3-Button Autoclavable Handpiece	Non-Sterile Sheath	Blunt & Sharp Dermal Tips
		
A901	A910 (10 Included)	A806 & A804 (50 each)

Disposable Electrodes	Solid Disposable Grounding Pad	Reusable Grounding Cord	Wall Mount Kit
			
ES01 (3 Included), ES02, ES20	A1202 (5 Included)	A1202C	A837

OTHER OPTIONAL ACCESSORIES

Footswitch	Reusable Bipolar Forceps*	Forcep Cord
		
A1203	A825 & A823	A827EU

Disposable and Reusable Electrodes*			Sterile Sheath
			
Lletz Loop Electrodes*	Resistick II™ Coated Electrodes*	Tungsten Loops*	A910ST

*More styles and accessories available — contact your Bovie representative for details.
Industry-standard accessories fit most electrosurgical generators

TECHNICAL SPECIFICATIONS

Technical Specifications...

Dimensions:		Line Voltage:	Line Frequency:
Height	16.0 cm (6.3 in.)	A950: 120 VAC A950-220: 240 VAC	50 - 60 Hz
Width	22.8 cm (9 in.)		
Depth	11.4 cm (4.5 in.)		
Weight	<3.2 kg (< 7 lbs)		
Monopolar	Output Power	Output Frequency	Rep Rate
Cut	60 W @ 500 Ω	357 kHz ± 50 kHz	N/A
Blend	60 W @ 800 Ω	357 kHz ± 50 kHz	30 kHz ± 5 kHz
Coagulation	60 W @ 1000 Ω	375 kHz ± 50 kHz	60 kHz ± 5 kHz
Fulguration	35 W @ 1000 Ω	575 kHz ± 50 kHz	30 kHz ± 5 kHz
Bipolar	30 W @ 200 Ω	520 kHz ± 14 kHz + 50 kHz	19 kHz ± 5 kHz
Warranty: limited two (2) year on unit only			

A950 Unit

A808T Tray
(Optional)

A812 Stand
(Optional)



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie® and Aaron® are registered trademarks of Bovie Medical Corporation